



WaterSaver

INNOVATIVE PRODUCTS FOR RESEARCH




NEL *direct*
online industrial & laboratory furniture
888.635.2949 | NELdirect.com

WaterSaver Faucet
701 W Erie St
Chicago, IL 60654

312 666 5500 TELEPHONE
312 666 5501 FACSIMILE
wsflab.com

L4160CV-158WSA

LABORATORY LIFT/TURN BALL VALVE ASSEMBLY, PANEL MOUNTED SINGLE, CHECK VALVE



NOTES:



1. LABORATORY BALL VALVE ASSEMBLY IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ANSI Z21.15 AND CGA 9.1.
2. HANDLE LOCKS IN CLOSED POSITION. LIFT HANDLE TO UNLOCK AND OPEN VALVE.
3. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED AT 150 PSI AIR PRESSURE. MAXIMUM WORKING PRESSURE IS 100 PSI.
4. SPECIFY IF CLEANING FOR OXYGEN AND HIGH PURITY GASES IS REQUIRED
5. FIXTURE HAS A POLISHED CHROME PLATED FINISH.

MEASUREMENTS MAY VARY $\pm 1/4"$ (6mm).

Drawing Number: **NELRS-12**

Revision Number: **090619-TEC**